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

CANDIDA AURIS

In Vitro Bactericidal Effects of 625, 525, and 425 nm Wavelength (Red, Green, and Blue) Light-Emitting Diode Irradiation

SangWoo Kim, MS,^{1*} JiSun Kim, PhD,^{1*} WonBong Lim, PhD,¹ SangMi Jeon, BS,¹ OkSu Kim, DDS, PhD,² Jeong-Tae Koh, DDS, PhD,³ Chang-Su Kim, PhD,⁴ HongRan Choi, DDS, PhD,¹ and OkJoon Kim, DDS, PhD¹

Abstract

Objective: The purpose of this study was to evaluate the relationship of 625, 525, and 425 nm wavelengths, providing average power output and effects on three common pathogenic bacteria. **Background data:** Ultraviolet (UV) light kills bacteria, but the bactericidal effects of UV may not be unique, as 425 nm produces a similar effect. The bactericidal effects of light-emitting diode (LED) wavelengths such as 625 and 525 nm have not been described. Before conducting clinical trials, the appropriate wavelength with reasonable dose and exposure time should be established. **Materials and methods:** The bactericidal effects of 625, 525, and 425 nm wavelength LED irradiation were investigated *in vitro* for the anaerobic bacterium *Porphyromonas gingivalis* and two aerobes (*Staphylococcus aureus* and *Escherichia coli* DH5α). Average power output was 6 mW/cm² for 1 h. The bacteria were exposed to LED irradiation for 1, 2, 4, and 8 h (21.6, 43.2, 86.4, and 172.8 J/cm², respectively). LED irradiation was performed during growth on agar and in broth. Control bacteria were incubated without LED irradiation. Bacterial growth was expressed in colony-forming units (CFU) and at an optical density at 600 nm in agar and broth. **Results:** The bactericidal effect of LED phototherapy depended upon wavelength, power density, bacterial viable number, and bacteria species. The bactericidal effect of 425 and 525 nm irradiation varied depending upon the bacterial inoculation, compared with unirradiated samples and samples irradiated with red light. Especially, *P. gingivalis* and *E. coli* DH5α were killed by 425 nm, and *S. aureus* growth was inhibited by 525 nm. However, the wavelength of 625 nm was not bactericidal for *P. gingivalis*, *E. coli* DH5α, or *S. aureus*. **Conclusions:** Irradiation at 625 nm light was not bactericidal to *S. aureus*, *E. coli*, and *P. gingivalis*, whereas wavelengths of 425 and 525 nm had bactericidal effects. *S. aureus* was also killed at 525 nm.

In Vitro Bactericidal Effects of 405-nm and 470-nm Blue Light

J. STEPHEN GUFFEY, Ed.D.,¹ and JAY WILBORN, M.Ed.²

ABSTRACT

Objective: The aim of this study was to determine the bactericidal effect of 405- and 470-nm light on two bacteria, *Staphylococcus aureus* and *Pseudomonas aeruginosa*, *in vitro*. **Background Data:** It is well-known that UV light kills bacteria, but the bactericidal effects of UV may not be unique since recent studies indicate that blue light produces a somewhat similar effect. The effects of blue light seem varied depending on wavelength, dose and the nature of the bacteria, hence this study. **Methods:** Two common aerobes, *Staphylococcus aureus* and *Pseudomonas aeruginosa*, and anaerobic *Propionibacterium acnes* were tested. Each organism was treated with Super Luminous Diode probes with peak emission at 405 and 470 nm. Treatment was timed to yield 1, 3, 5, 10, and 15 J cm⁻² doses. Colony counts were performed and compared to untreated controls. **Results:** The 405-nm light produced a dose dependent bactericidal effect on *Pseudomonas aeruginosa* and *Staphylococcus aureus* ($p < .05$), achieving as much as 95.1% and nearly 90% kill rate for each, respectively. The 470-nm light effectively killed *Pseudomonas aeruginosa* at all dose levels, but only killed *Staphylococcus aureus* at 10 and 15 J cm⁻². With this wavelength, as much as 96.5% and 62% reduction of *Pseudomonas aeruginosa* and *Staphylococcus aureus* was achieved, respectively. Neither of the two wavelengths proved bactericidal with anaerobic *Propionibacterium acnes*. **Conclusion:** The results indicate that, *in vitro*, 405- and 470-nm blue light produce dose dependent bactericidal effects on *Pseudomonas aeruginosa* and *Staphylococcus aureus* but not *Propionibacterium acnes*.

J. Periodontol. Implant Sci. 2003 Apr;32(5):110-119
http://dx.doi.org/10.1007/s10000000046
pISSN 1003-0279 eISSN 1003-0596

Research Article Periodontal Science

In vitro investigation of the antibacterial and anti-inflammatory effects of LED irradiation

Jungwon Lee , Hyun-Yong Song , Sun-Hee Ahn , Woosub Song , Yang-Jo Seol , Yong-Moo Lee , Ki-Tae Koo

¹Department of Periodontology, School of Dentistry and Dental Research Institute, Seoul National University, Seoul, Korea
²Oral Stop Specialty Center, Seoul National University Dental Hospital, Seoul, Korea
³Medical & Bio Photonics Research Center, Korea Photonics Technology Institute (KOPTI), Gwangju, Korea

OPEN ACCESS

Received: Feb 21, 2020
Revised: Mar 16, 2020
Accepted: Apr 4, 2020
Published online: Sep 1, 2020

*Correspondence:

Yong-Moo Lee
Department of Periodontology and Dental Research Institute, School of Dentistry, Seoul National University, 101 Daehak-ro, Jongno-gu, Seoul 05150, Korea
Email: yml@plaza.snu.ac.kr
Tel: +82-2-762-3034
Fax: +82-2-764-0051

Ki-Tae Koo

Department of Periodontology and Dental Research Institute, School of Dentistry, Seoul National University, 101 Daehak-ro, Jongno-gu, Seoul 05150, Korea
Email: ktm@plaza.snu.ac.kr
Tel: +82-2-762-3038
Fax: +82-2-764-0051

Copyright © 2020 Korean Academy of Periodontology
This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<https://creativecommons.org/licenses/by-nc/4.0/>)

ABSTRACT

Purpose: This study aimed to investigate the proper wavelengths for safe levels of light-emitting diode (LED) irradiation with bactericidal and photomodulation effects *in vitro*. **Methods:** Cell viability tests of fibroblasts and osteoblasts after LED irradiation at 470, 525, 590, 630, and 850 nm were performed using the thiazolyl blue tetrazolium bromide assay. The bactericidal effect of 470-nm LED irradiation was analyzed with *Streptococcus pyogenes*, *Aggregatibacter actinomycetemcomitans*, *Pasteurella acnes*, *Porphyromonas gingivalis*, and *Neisseria meningitidis*. Levels of nitric oxide, a proinflammatory mediator, were measured to identify the anti-inflammatory effect of LED irradiation on lipopolysaccharide-stimulated inflammation in RAW 264.7 macrophages. **Results:** LED irradiation at wavelengths of 470, 525, 590, 630, and 850 nm showed no cytotoxic effect on fibroblasts and osteoblasts. LED irradiation at 630 and 850 nm led to fibroblast proliferation compared to no LED irradiation. LED irradiation at 470 nm resulted in bactericidal effects on *S. pyogenes*, *A. actinomycetemcomitans*, *P. acnes*, *P. gingivalis*, and *N. meningitidis*. Lipopolysaccharide (LPS)-induced RAW 264.7 inflammation was reduced by irradiation with 525-nm LED before LPS treatment and irradiation with 630-nm LED after LPS treatment; however, the effects were limited. **Conclusions:** LED irradiation at 470 nm showed bactericidal effects, while LED irradiation at 525 and 630 nm showed preventive and treatment effects on LPS-induced RAW 264.7 inflammation. The application of LED irradiation has potential as an adjuvant in periodontal therapy, although further investigations should be performed *in vivo*.

Keywords: Bacteria; Fibroblasts; Inflammation; Osteoblasts; Photobiology; Phototherapy

scientific reports

OPEN

Efficacy of violet-blue light to inactive microbial growth

Davide Amodeo¹, Valentina Lucarelli², Isa De Palma¹, Alessandro Puccio¹, Nicola Nante³, Gabriele Cevenini¹ & Gabriele Messina^{2,3}

The increase in health care-associated infections and antibiotic resistance has led to a growing interest in the search for innovative technologies to solve these problems. In recent years, the interest of the scientific community has focused on violet-blue light at 405 nm (VBL405). This study aimed to assess the VBL405 efficiency in reducing microbial growth on surfaces and air. This descriptive study run between July and October 2020. Petri dishes were contaminated with *P. aeruginosa*, *E. coli*, *S. aureus*, *S. typhimurium*, *K. pneumoniae* and were placed at 2 and 3 m from a LED light source having a wavelength peak at 405 nm and an irradiance respectively of 967 and 497 μW/cm². Simultaneously, the air in the room was sampled for 5 days with two air samplers (SAS) before and after the exposure to the VBL405 source. The highest microbial reduction was reached 2 m directly under the light source: *S. typhimurium* (2.93 log₁₀), *K. pneumoniae* (2.30 log₁₀), *S. aureus* (3.98 log₁₀), *E. coli* (3.83 log₁₀), *P. aeruginosa* (3.86 log₁₀). At a distance of 3 m from the light source, the greatest reduction was observed for *S. aureus* (3.49 log₁₀), and *P. aeruginosa* (3.80 log₁₀). An average percent microbial reduction of about 70% was found in the sampled air after 12 h of exposure to VBL405. VBL405 has proven to contrast microbial growth on the plates. Implementing this technology in the environment to provide continuous disinfection and to control microbial presence, even in the presence of people, may be an innovative solution.

HHS Public Access

Author manuscript
Lasers Surg Med. Author manuscript; available in PMC 2016 April 16.

Published in final edited form as:
Lasers Surg Med. 2015 March; 47(3): 266–272. doi:10.1002/lsm.22327.

Article

Determination of Photosensitizing Potential of Lapachol for Photodynamic Inactivation of Bacteria

Regiane G. Lima , Raphael S. Flores, Gabriella Miessi, Jhoenne H. V. Pulcherio, Laís F. Aguilera, Leandro O. Araujo, Samuel L. Oliveira and Anderson R. L. Calves

Optics and Photonics Group, Institute of Physics, Federal University of Mato Grosso do Sul, P.O. Box 549, 79070-900 Campo Grande, MS, Brazil; regigodoy@gmail.com (R.G.L.); raphael.rflores@gmail.com (R.S.F.); gabriella.miessi@ufms.br (G.M.); jhoenne_vasconcelos@ufms.br (J.H.V.P.); lais_fernandes13@hotmail.com (L.F.A.); leandro.oliveira22@outlook.com (L.O.A.); samuel.oliveira@ufms.br (S.L.O.)

* Correspondence: anderson.calves@ufms.br

Abstract: Antimicrobial photodynamic inactivation (aPDI) offers a promising alternative to combat drug-resistant bacteria. This study explores the potential of lapachol, a natural naphthoquinone derived from *Tabeaia avellanae*, as a photosensitizer (PS) for aPDI. Lapachol's photosensitizing properties were evaluated using *Staphylococcus aureus* and *Escherichia coli* strains under blue LED light (450 nm). UV-vis spectroscopy confirmed lapachol's absorption peak at 482 nm, aligning with effective excitation wavelengths for phototherapy. Photoinactivation assays demonstrated significant bacterial growth inhibition, achieving complete eradication of *S. aureus* at 25 μg mL⁻¹ under light exposure. Scanning electron microscopy (SEM) revealed morphological damage in irradiated bacterial cells, confirming lapachol's bactericidal effect. This research underscores lapachol's potential as a novel photosensitizer in antimicrobial photodynamic therapy, addressing a critical need in combating antibiotic resistance.

Optimization of the Antimicrobial Effect of Blue Light on Methicillin-Resistant *Staphylococcus aureus* (MRSA) *In Vitro*

Violet V. Bumah, PhD¹, Daniela S. Masson-Meyers, PhD¹, Susan Cashion, PhD¹, and Chukuka S. Enwemeka, PhD^{1,2,*}

¹College of Health Sciences, University of Wisconsin-Milwaukee, Milwaukee, Wisconsin 53211
²San Diego State University, San Diego, California 92182

Abstract

Background and Objective—In previous studies, we showed that irradiation with 405 nm or 470 nm light suppresses up to 92% methicillin-resistant *Staphylococcus aureus* (MRSA) growth *in vitro* and that the remaining bacteria re-colonize. In this study, the aim was to develop a protocol that yields 100% MRSA growth suppression.

Materials and Methods—We cultured 3×10^6 and 5×10^6 CFU/ml USA300 strain of MRSA and then irradiated each plate with varying fluences of 1–60 J/cm² of 405 nm or 470 nm light, either once or twice at 6 hours intervals. Next, we plated 7×10^6 CFU/ml and irradiated it with 45, 50, 55, or 60 J/cm² fluence, once, twice, or three at the same 6 hours intervals. In a third experiment, the same culture density was irradiated with 0, 165, 180, 220, or 240 J/cm², either once, twice, or three.

Results—Irradiation with either wavelength significantly reduced the bacterial colonies regardless of bacterial density ($P < 0.05$). At 3×10^6 CFU/ml density, nearly 40% and 50% growth of MRSA were suppressed with as little as 3 J/cm² of 405 nm and 470 nm wavelengths, respectively. Moreover, 100% of the colonies were suppressed with a single exposure to 55 or 60 J/cm² of 470 nm light or double treatment with 50, 55, or 60 J/cm² of 405 nm wavelength. At 5×10^6 CFU/ml density, irradiating twice with 50, 55, or 60 J/cm² of either wavelength suppressed bacterial growth completely, lower fluences did not. The denser 7×10^6 CFU/ml culture required higher doses to achieve 100% suppression, either one shot with 220 J/cm² of 470 nm light or two shots of the same dose using 405 nm.

Conclusion—The bactericidal effect of blue light can be optimized to yield 100% bacterial growth suppression, but with relatively high fluences for dense bacterial cultures, such as 7×10^6 CFU/ml.